



Model 6530 Series Racks for Gas Cylinder Storage



Description

The basic models secure from 2 to 6 full size cylinders (up to 10" diameter) per rack. These racks are constructed of welded steel and can be anchored in concrete with 1/2" steel anchor bolts*. Models come with double restraints. All racks are finished in a tough exterior grade baked polyurethane finish.

* Not supplied.

Ordering Information

Model 6530 Series

Model No.	No. of Cylinders	No. of Restraints	Rack Height	Rack Width	Rack Depth	Shipping Weight
6531-2	2	2	30"	30"	16"	91 lbs
6532-2	4	2	30"	30"	26"	101 lbs
6533-2	6	2	30"	30"	38"	136 lbs