




## Model 700, 701, and 702 Small Cylinder Stand



### Description

#### Model 700

Small Cylinder Stand. A safe and economical product for use with 4 inch diameter gas bottles commonly used for medical oxygen and in laboratories. The stand helps prevent the gas bottles from tipping over, locking the gas bottles in place with nylon set screws.

#### Specifications

Size: 9 x 9 x 7  
 Shipping Weight: 3 lbs

#### Model 701

Ring stand for medium 5 to 6.75 inch gas bottles. Four nylon set screws.

#### Specifications

Size: 15 x 15 x 7  
 Shipping Weight: 5 lbs

#### Model 702

Ring stand for medium 7 to 8 inch gas bottles (3I and 3R), and maximum cylinder height of 20 inches. Four nylon set screws.

#### Specifications

Size: 17 x 17 x 7  
 Shipping Weight: 8 lbs

### Ordering Information

Model No.	Description
700	Steel Stand for 3.25 to 4.25" Bottles
701	Ring Stand for 5 to 6.75" Bottles
702	Ring Stand for 7 to 8" Bottles (3I and 3R) 20" Max Cylinder Height
HOL-0026-XX	Mini-MAT 6R Stand

### Other Cylinder Holder Equipment

**Lecture Bottle Holders:** See page 433 'Lecture Bottle Equipment' for product information.