



PUR-Gas™ In-Line Purifier Systems (continued)

Install and Replace Traps without Tools and in Seconds!



1 Align End Connector into Purifier



2 Screw to Tighten

Ordering Information

STANDARD IN-LINE "TRAP" Models

Model No.	Description
PUR-IL-MT1	PUR-Gas In-Line Moisture (H ₂ O) Trap
PUR-IL-OT1	PUR-Gas In-Line Oxygen (O ₂) Trap
PUR-IL-HT1	PUR-Gas In-Line Hydrocarbon (HC) Trap
PUR-IL-OMT1	PUR-Gas In-Line Combi Oxygen/Moisture Trap
PUR-IL-TRT1	PUR-Gas In-Line Triple Oxygen/Moisture/HC Trap
PUR-IL-ITRT1	PUR-Gas "Breakthrough" In-Line (He filled) Triple Oxygen/Moisture/HC Indicator*
PUR-IL-TRT2	PUR-Gas In-Line Triple Oxygen/Moisture/HC Trap (He filled)

* "Indicating" Traps provide visual indication for determining the optimum time for replacement and maximizing operating effectiveness. (a) The housing for the "Indicating" In-Line Triple Trap is made from glass encased with a polycarbonate coating; all other in-line trap housings are stainless steel material and are "Non-Indicating." (b) The term "He filled" applies to those specific traps that have been purged and filled internally with helium gas; all remaining trap models are filled with argon gas and should be purged after installation.

"END-to-END" Tube Connectors ("1" Set required for each In-Line Purifier purchased)

Model No.	Description
PUR-IL-CBR2	1/8" Brass In-Line Connector (set of 2)
PUR-IL-CBR4	1/4" Brass In-Line Connector (set of 2)
PUR-IL-CSS2	1/8" Stainless Steel In-Line Connector (set of 2)
PUR-IL-CSS4	1/4" Stainless Steel In-Line Connector (set of 2)
PUR-IL-DTCSS4	1/4" Stainless Steel SERIES In-Line Connector (qty 1)**

**The "SERIES" In-Line Connector is a single connector ONLY used for connecting '2' or more traps in a direct END-to-END series installation arrangement where it is desired to have gas flow from one trap immediately (sequentially) into a second (or even third) trap connected in "SERIES" with the first trap.

OPTIONAL & REPLACEMENT Parts

Model No.	Description
PUR-IL-WMC4	Wall-Mount Clamp Set (for Standard Models)
PUR-IL-ORS10	Spare O-Ring Set for In-Line Connectors (pkg of 10)

PRODUCT APPLICATION RECOMMENDATIONS:

1. END-to-END Connectors are primarily intended for use with stainless steel and brass tubing; connections to plastic tubing can be made; however, it is recommended that a metal insert be utilized at each inlet/outlet connection point where plastic tubing connections are required.
2. The proper installation orientation for all purifiers is in the "vertical" position.
3. Mounting screws for Wall-Mount Clamps are not included.



In-line Purifier with Stainless Steel Housing



1/8" or 1/4" Fittings on End Connectors



In-Line Triple Purifier with Visual Indication



In-Line Purifier with Wall-Mount Clips