



Model 500 Small Cylinder Stand



Description

This Cylinder Stand supports the number 3 and 4 size cylinders. It is constructed of nickel plated steel and is a simple, safe, economical device to stabilize the position of small cylinders. It accommodates cylinders with diameters from 4" to 7-3/8".

Specifications

Size: 12-3/4"W x 11-1/8"H
Shipping Weight: 4 lbs

Ordering Information

Model No.	Description
500	Small Cylinder Stand
HOL-0026-XX	Mini-MAT 6R Stand

Other Cylinder Holder Equipment

Lecture Bottle Holders: See page 433 'Lecture Bottle Equipment' for product information.