


## Models 504 & 505 Lecture Bottle Holders



Model 505



Model 504

### Description

Lecture bottles have rounded ends and require some means of support when in use. The Model 505 Non-Tip Stand is recommended when a regulator or other equipment is attached to the lecture bottle. The Model 504 Wall Mount Bracket is made of plastic coated steel and has spring steel clamps to ensure firm, safe support of the lecture bottle.

### Ordering Information

Model No.	Type	Shipping Weight
504	Wall Mount Bracket	2 lbs
505	Non-Tip Stand (Aluminum-Satin Finish)	2 lbs

### Other Lecture Bottle Equipment

**Lecture Bottle Control Valves:** See page 408 'Valves' for product information.

## Model 1191 Lecture Bottle Storage Cabinet



### Description

The Model 1191 Lecture Bottle Storage Cabinet ensures safe storage of hazardous gases in lecture bottles. Two (2) lecture bottle holders keep the lecture bottles safely in place. Minimum air velocity of 77 ft/min at 100 SCFM ensures proper venting to prevent any possible leaks from reaching the work environment. The compact size of the Cabinet allows for placement at the lab bench for easy access and handling. Safety glass window, key lock and NEMA 4 enclosure rating ensure additional safety. Equipped with two (2) Model 504 Cylinder Holders. Not designed for Lecture Bottles in service.

### Construction:

Body:	14 Gauge Cold Rolled Steel
Door:	14 Gauge Cold Rolled Steel
Seams:	Continuously welded
Gasket:	Oil resistant
Finish:	Grey polyester powder
Internal Finish:	White enamel
Lecture Bottle Holders:	2
(Model 504)	
Shipping Weight:	50 lbs

### Dimensions:

Cabinet - External	
Height:	24"
Width:	20"
Depth:	8-7/8"
Stack	
Diameter:	4"
Height:	2"

### Ordering Information

Model No.	Description
1191	Lecture Bottle Storage Cabinet